

SINCE 1978



OPTICS TECHNOLOGY

Manufacturer & Exporters of : LABORATORY TESTING INSTRUMENTS



An ISO 9001:2008 CERTIFIED
An ISO 10002:2004 CERTIFIED

CE Certified

ABBE Refractometer with Imported Optics



ABBE Refractometer with Imported Optics

Clarity is more compared to normal ABBE refractometer.
Background of Scale is green with clear scale, divisions in black.
Both x scale & measuring scale in single eye piece.

Measuring range: nD 1.3000 to 1.7000
Sugar : 0 to 95% with both refractive index Sugar Scale

2/3, DSIDC, Community Works Centre, J-Block, Shakurpur, (Opp. Britannia Ind. Ltd.) Delhi-110034 (India)

Phone: 0091-11-27101383, 27184490, Fax: 0091-11-27108890

E-mail: info@opticstechnology.in, opticstechnology@gmail.com

Website: www.opticstechnology.in