

SINCE 1978



# OPTICS TECHNOLOGY

Manufacturer & Exporters of : LABORATORY TESTING INSTRUMENTS



An ISO 9001:2008 CERTIFIED  
An ISO 10002:2004 CERTIFIED

**CE Certified**

## ABBE Refractometer



ABBE Refractometer

ABBE refractometer is a standard instrument for measuring refractive index and dispersion of liquids and solids and also for determining the percentage contents of sugar in water solutions.

Prism box is water jacketed for circulation of water at constant temperature if required, for precise measurement.

**Measuring Brix range :** Refractive Index : nD 1.3000 to 1.7000

**Sugar :** 0 to 95% with both refractive index Sugar Scale.

2/3, DSIDC, Community Works Centre, J-Block, Shakurpur, (Opp. Britannia Ind. Ltd.) Delhi-110034 (India)

Phone: 0091-11-27101383, 27184490, Fax: 0091-11-27108890

E-mail: [info@opticstechnology.in](mailto:info@opticstechnology.in), [opticstechnology@gmail.com](mailto:opticstechnology@gmail.com)

Website: [www.opticstechnology.in](http://www.opticstechnology.in)