

SINCE 1978



OPTICS TECHNOLOGY

Manufacturer & Exporters of : LABORATORY TESTING INSTRUMENTS



An ISO 9001:2008 CERTIFIED
An ISO 10002:2004 CERTIFIED

CE Certified

Digital ABBE Refractometer



Digital ABBE Refractometer

Digital display of refractive index

Measurement Range: refractive index nd : 1.3000-1.7000

Measurement Range: Brix (Brix) : 0-95%

Accuracy: refractive index nd : $\leq \pm 0.0002$

Accuracy: Brix (Brix) : $\leq \pm 0.1\%$

Showing temperature range (minimum 0.1 °C) : 0-50 °C

2/3, DSIDC, Community Works Centre, J-Block, Shakurpur, (Opp. Britannia Ind. Ltd.) Delhi-110034 (India)

Phone: 0091-11-27101383, 27184490, Fax: 0091-11-27108890

E-mail: info@opticstechnology.in, opticstechnology@gmail.com

Website: www.opticstechnology.in