

SINCE 1978



OPTICS TECHNOLOGY

Manufacturer & Exporters of : LABORATORY TESTING INSTRUMENTS



An ISO 9001:2008 CERTIFIED
An ISO 10002:2004 CERTIFIED

CE Certified

Automatic Polarimeter



Automatic Polarimeter

TECHNICAL SPECIFICATIONS:

Measurement Range	45° ~ +45°
Accuracy	±0.02°
Light Source	Sodium lamp
Optical Wavelengths	589.44nm
Stable Time	5 minutes
Tube Length	up to 200mm
Power Requirements	220VAC/50Hz

Automatic polarimeter is the basic scientific instrument used to measure angle of rotation, specific optical rotation of optically active substance, built in auto controlling system and displays the optical rotation directly within a short duration as soon as the sample filled tube is placed.

2/3, DSIDC, Community Works Centre, J-Block, Shakurpur, (Opp. Britannia Ind. Ltd.) Delhi-110034 (India)

Phone: 0091-11-27101383, 27184490, Fax: 0091-11-27108890

E-mail: info@opticstechnology.in, opticstechnology@gmail.com

Website: www.opticstechnology.in