

SINCE 1978



# OPTICS TECHNOLOGY

Manufacturer & Exporters of : LABORATORY TESTING INSTRUMENTS



An ISO 9001:2008 CERTIFIED  
An ISO 10002:2004 CERTIFIED

**CE Certified**

## Bead Sterilizer



Bead Sterilizer

Glass Bead sterilizers are very convenient to use over alcohol / flame for sterilizing implements (e.g. forceps, scalpels etc..) in the laboratory.

### FEATURES:

Quick and clean alternative to conventional methods of sterilization, ensure total sterilization by destruction of anti micro organisms.

Operating temperature upto 250-265°C.

Fully insulated body.

The unit is compact and can be placed on any Laminar Air Flow Work tables in a clean room atmosphere.

2/3, DSIDC, Community Works Centre, J-Block, Shakurpur, (Opp. Britannia Ind. Ltd.) Delhi-110034 (India)

Phone: 0091-11-27101383, 27184490, Fax: 0091-11-27108890

E-mail: [info@opticstechnology.in](mailto:info@opticstechnology.in), [opticstechnology@gmail.com](mailto:opticstechnology@gmail.com)

Website: [www.opticstechnology.in](http://www.opticstechnology.in)