

SINCE 1978



OPTICS TECHNOLOGY

Manufacturer & Exporters of : LABORATORY TESTING INSTRUMENTS



An ISO 9001:2008 CERTIFIED
An ISO 10002:2004 CERTIFIED

CE Certified

Bulk Density Apparatus



Bulk Density Apparatus

Digital Automatic Tap Density Test Apparatus is designed to meet the requirements of USP Standards.

It can hold two cylinders.

The mechanism provides simultaneously rotating and tapping motion to the cylinders for packing the sample under test evenly and uniformly together.

Body is made of M.S., duly powder coated / S.S. body (as per GMP requirement).

The user has to select the No. of stokes in the counter & after selecting, just press start switch, it will start the stokes & the Digital display shows the No. of stokes.

Cylinder Capacity: 100 ml

2/3, DSIDC, Community Works Centre, J-Block, Shakurpur, (Opp. Britannia Ind. Ltd.) Delhi-110034 (India)

Phone: 0091-11-27101383, 27184490, Fax: 0091-11-27108890

E-mail: info@opticstechnology.in, opticstechnology@gmail.com

Website: www.opticstechnology.in