

SINCE 1978



OPTICS TECHNOLOGY

Manufacturer & Exporters of : LABORATORY TESTING INSTRUMENTS



An ISO 9001:2008 CERTIFIED
An ISO 10002:2004 CERTIFIED

CE Certified

Butyro Refractometer



Butyro Refractometer

SUGAR AND OIL BUTYRO REFRACTOMETER)

This instrument is used to check the B.R. VALUE, in Oil.

Divisions on the oil scale are sub-divided by Micro meter screw

Error due to temperature changes in sugar solutions can be eliminated by temperature corrections device provided.

It is an indispensable tool for use of all industries using sugar, oil & many other similar substances.

Butyro Refractometer Index of Measuring range: 1.42 to 1.49 to an accuracy of 0.0002. Micrometer allows accuracy of 0.01% on the oil (Butyro) scale.

2/3, DSIDC, Community Works Centre, J-Block, Shakurpur, (Opp. Britannia Ind. Ltd.) Delhi-110034 (India)

Phone: 0091-11-27101383, 27184490, Fax: 0091-11-27108890

E-mail: info@opticstechnology.in, opticstechnology@gmail.com

Website: www.opticstechnology.in