

SINCE 1978



# OPTICS TECHNOLOGY

Manufacturer & Exporters of : LABORATORY TESTING INSTRUMENTS



An ISO 9001:2008 CERTIFIED  
An ISO 10002:2004 CERTIFIED

**CE Certified**

## Carbon Sulphur Apparatus



Carbon Sulphur Apparatus

Heavy duty polished wooden frame with best quality clamps to hold glass parts. Furnace outer body M.S. sheet duly powder coated having temp. control box with pyrometer / digital temp. controller. Supplied with glass parts absorption vessel & Acid Jar, Combustion Boat etc.

Available in two models:

- (i) 15 A control unit series
- (ii) 35 A control unit parallel

2/3, DSIDC, Community Works Centre, J-Block, Shakurpur, (Opp. Britannia Ind. Ltd.) Delhi-110034 (India)

Phone: 0091-11-27101383, 27184490, Fax: 0091-11-27108890

E-mail: [info@opticstechnology.in](mailto:info@opticstechnology.in), [opticstechnology@gmail.com](mailto:opticstechnology@gmail.com)

Website: [www.opticstechnology.in](http://www.opticstechnology.in)