

SINCE 1978



# OPTICS TECHNOLOGY

Manufacturer & Exporters of : LABORATORY TESTING INSTRUMENTS



An ISO 9001:2008 CERTIFIED  
An ISO 10002:2004 CERTIFIED

**CE Certified**

## Chemical Storage Cabinet



Chemical Storage Cabinet

These specially designed cabinets provide a wash of sterile filtered, clean air through HEPA Filters having an efficiency of 99.97% down to 0.3 micron particle size. This cabinet protects dust and other particulates on Bottles and chemicals. For use in the Laboratory for chemical storage.

Self standing cabinet provides down flow clean Air in the cabinet.

### FEATURES

Large interior capacity for handling various sizes Bottles / Chemicals. Wire Mesh Trays are Provided in the Chamber so that all the air will be fully recirculated. Low power consumption suitable for operation on single phase supply.

**Inner coated with Acid Resistant & Fire Retardent Resin and Outer M.S. Body Duly Powder Coated**

### Working Area Size (In Ft.)

2' x 1.5" x 6'  
3' x 1.5" x 6'  
4' x 1.5" x 6'

2/3, DSIDC, Community Works Centre, J-Block, Shakurpur, (Opp. Britannia Ind. Ltd.) Delhi-110034 (India)

Phone: 0091-11-27101383, 27184490, Fax: 0091-11-27108890

E-mail: [info@opticstechnology.in](mailto:info@opticstechnology.in), [opticstechnology@gmail.com](mailto:opticstechnology@gmail.com)

Website: [www.opticstechnology.in](http://www.opticstechnology.in)