

SINCE 1978



OPTICS TECHNOLOGY

Manufacturer & Exporters of : LABORATORY TESTING INSTRUMENTS



An ISO 9001:2008 CERTIFIED
An ISO 10002:2004 CERTIFIED

CE Certified

COD Digestor



COD Digestor

COD Digestion Apparatus is a solid block heated unit with provision for 6 samples at a time in a Reaction Vessel. All the sample are heated at 150°C within + 1% accuracy with utmost uniformity of temperature.

These are fitted with digital temperature controllers, buzzer and timer. Further, our range is available in the following specifications:

- Heating block type
- Temperature range from 5 °C above and ambient to 150 °C \pm 1 °C
- Designed with glass reactions vessels and air condensers
- Supplied with stand to accommodate vessels
- Operate on 230 volts, 50 Hz, single phase
- Available sizes for 6 samples, 15 samples, 20 samples & 25 Samples at a time.