

SINCE 1978



# OPTICS TECHNOLOGY

Manufacturer & Exporters of : LABORATORY TESTING INSTRUMENTS



An ISO 9001:2008 CERTIFIED  
An ISO 10002:2004 CERTIFIED

**CE Certified**

## Digital Conductivity Meter



Digital Conductivity Meter

### Table Model With The Following Specifications :

Range : 0-200 microsemens , 0-2 milisemens , 0-20 milisemens, 0-200 milisemens, 0-1000 milisemens

Accuracy: +/- 1% FS +/-1

Digit Resolution: 0.01% microsemens per cm

Cell : 0.05 To 10.5 Adjustable

Provision of Temperature Compensation

Digital Cell Constant Adjustment Along With Conductivity Cell And Temperature Probe

Display : 3 ½ Digit LED

Power Supply: 230V +/- 10%, 50Hz +/-6%

2/3, DSIDC, Community Works Centre, J-Block, Shakurpur, (Opp. Britannia Ind. Ltd.) Delhi-110034 (India)

Phone: 0091-11-27101383, 27184490, Fax: 0091-11-27108890

E-mail: [info@opticstechnology.in](mailto:info@opticstechnology.in), [opticstechnology@gmail.com](mailto:opticstechnology@gmail.com)

Website: [www.opticstechnology.in](http://www.opticstechnology.in)