

SINCE 1978



OPTICS TECHNOLOGY

Manufacturer & Exporters of : LABORATORY TESTING INSTRUMENTS



An ISO 9001:2008 CERTIFIED
An ISO 10002:2004 CERTIFIED

CE Certified

Digital Butyro Refractometer



Digital Butyro Refractometer

No need of circulating water bath, displays result of ATC value only

Features:-

- * Calibration with Standard liquid
- * Large, easy to read display
- * External Light Interference (ELI)
- * Display of both Butyro & RI scale
- * Temperature Compensated scale at 40° C
- * Automatic & Accurate (Accuracy : Butyro ± 0.5)

Measurement Range	Butyro 30.0 to 90.0
(Converted at 40°C)	RI 1.4450 to 1.4850
Minimum Indication	Butyro 0.1 RI 0.0001
Measurement Accuracy	Butyro ± 0.5 (at 40°C) RI ± 0.0003 (at 40°C)
Measurement Temperature	10 to 50°C (Automatic Temperature Compensation)

2/3, DSIDC, Community Works Centre, J-Block, Shakurpur, (Opp. Britannia Ind. Ltd.) Delhi-110034 (India)

Phone: 0091-11-27101383, 27184490, Fax: 0091-11-27108890

E-mail: info@opticstechnology.in, opticstechnology@gmail.com

Website: www.opticstechnology.in