

SINCE 1978



# OPTICS TECHNOLOGY

Manufacturer & Exporters of : LABORATORY TESTING INSTRUMENTS



An ISO 9001:2008 CERTIFIED  
An ISO 10002:2004 CERTIFIED

**CE Certified**

## Digital Conductivity Meter



Digital Conductivity Meter

Table Model With The Following Specifications :

Range : 0-200 microsemens , 0-2 milisemens , 0-20 milisemens, 0-200 milisemens, 0-1000 milisemens  
Accuracy: +/- 1% FS +/-1  
Digit Resolution: 0.01% microsemens per cm  
Cell : 0.05 To 10.5 Adjustable  
Provision of Temperature Compensation  
Digital Cell Constant Adjustment Along With Conductivity Cell And Temperature Probe  
Display : 3 ½ Digit LED  
Power Supply: 230V +/- 10%, 50Hz +/-6%