

SINCE 1978



# OPTICS TECHNOLOGY

Manufacturer & Exporters of : LABORATORY TESTING INSTRUMENTS



An ISO 9001:2008 CERTIFIED  
An ISO 10002:2004 CERTIFIED

**CE Certified**

## Digital pH Meter



Digital pH Meter

- Auto Temperature Compensation Facility
- 3 ½ digit LED Display
- Auto Polarity & Decimal Indication
- Slope Control

DISPLAY	3 ½ Digit LED	
RANGE	pH	0 to 14.00
	mV	0 to ± 1999
RESOLUTION	pH	0.01
	mV	± 1 digit
ACCURACY	pH	± 0.01 pH
	mV	± 1 digit
TEMPERATURE	Auto: 0 to 100oC	
COMPENSATION	Manual : 0 to 100oC	
SLOPE CONTROL	80 to 120%	
POWER	230V ± 10% AC, 50 Hz	

2/3, DSIDC, Community Works Centre, J-Block, Shakurpur, (Opp. Britannia Ind. Ltd.) Delhi-110034 (India)

Phone: 0091-11-27101383, 27184490, Fax: 0091-11-27108890

E-mail: [info@opticstechnology.in](mailto:info@opticstechnology.in), [opticstechnology@gmail.com](mailto:opticstechnology@gmail.com)

Website: [www.opticstechnology.in](http://www.opticstechnology.in)