

SINCE 1978



OPTICS TECHNOLOGY

Manufacturer & Exporters of : LABORATORY TESTING INSTRUMENTS



An ISO 9001:2008 CERTIFIED
An ISO 10002:2004 CERTIFIED

CE Certified

Digital Refractometer

Digital Indication of Brix and temperature



Digital Refractometer

(i)	Range	0 to 85% Brix	0 to 60°C (32 to 178°F)
(ii)	Resolution	± 0.1% Brix	± 0.1°C (0.1°F)
(iii)	Accuracy	± 0.2% Brix	± 0.3°C (0.5°F)
(iv)	Temperature Compensation	Automatic Between 10 and 40°C (50 - 104°F)	
(v)	Measurement Time	Approximately 15 seconds	

2/3, DSIDC, Community Works Centre, J-Block, Shakurpur, (Opp. Britannia Ind. Ltd.) Delhi-110034 (India)

Phone: 0091-11-27101383, 27184490, Fax: 0091-11-27108890

E-mail: info@opticstechnology.in, opticstechnology@gmail.com

Website: www.opticstechnology.in