

SINCE 1978



OPTICS TECHNOLOGY

Manufacturer & Exporters of : LABORATORY TESTING INSTRUMENTS



An ISO 9001:2008 CERTIFIED
An ISO 10002:2004 CERTIFIED

CE Certified

DIGITAL SPECTROPHOTOMETER



DIGITAL SPECTROPHOTOMETER

VISIBLE RANGE

Range	340 To 960 nm
Output	% T : 0 To 100 %, Abs: 0 To 1.999, Conc. : 0 To 1999
Resolution	0.1 % T, Abs :0.001 Conc. :1
Detection	Silicon Photodiode
Spectral Bandwidth	5 Nm
Wave Length Resolution	1 Nm
Wave Length Accuracy	± 2.5 Nm
Repeatability	± 2 Nm
Display	Data : 3 1/2 Digit Led
Wavelength	4 Digit LED Display

NOTE: MICROPROCESSOR BASED UV/VIS SPECTROPHOTOMETERS, UV SPECTROPHOTOMETERS & DOUBLE BEAM SPECTROPHOTOMETERS ARE ALSO AVAILABLE ON REQUEST

2/3, DSIDC, Community Works Centre, J-Block, Shakurpur, (Opp. Britannia Ind. Ltd.) Delhi-110034 (India)

Phone: 0091-11-27101383, 27184490, Fax: 0091-11-27108890

E-mail: info@opticstechnology.in, opticstechnology@gmail.com

Website: www.opticstechnology.in