

SINCE 1978



OPTICS TECHNOLOGY

Manufacturer & Exporters of : LABORATORY TESTING INSTRUMENTS



An ISO 9001:2008 CERTIFIED
An ISO 10002:2004 CERTIFIED

CE Certified

Digital TDS Meter



Digital TDS Meter

Table Model With The Following Specifications :

Range : 0-200 PPM , 0-2000 PPM, 0-20 PPT, 0-200 PPT, 0-1000 PPT
Accuracy: +/- 1% FS +/-1
Cell : 0.05 To 10.5 Adjustable
Provision of Manual Temperature Compensation
Digital Cell Constant Adjustment TDS cell
Display : 3 ½ Digit LED
Power Supply: 230V +/- 10%, 50Hz +/-6%

ALSO AVAILABLE MICROPROCESSOR BASED CONDUCTIVITY /
TDS METER

2/3, DSIDC, Community Works Centre, J-Block, Shakurpur, (Opp. Britannia Ind. Ltd.) Delhi-110034 (India)

Phone: 0091-11-27101383, 27184490, Fax: 0091-11-27108890

E-mail: info@opticstechnology.in, opticstechnology@gmail.com

Website: www.opticstechnology.in