

SINCE 1978



OPTICS TECHNOLOGY

Manufacturer & Exporters of : LABORATORY TESTING INSTRUMENTS



An ISO 9001:2008 CERTIFIED
An ISO 10002:2004 CERTIFIED

CE Certified

Front Loading Table Top Autoclave



- Front Loading Autoclave is fully automatic front loading autoclave with dry cycle.
- Electronic Circuit inside the autoclave with programmable integrated circuit chip.
- With dry cycle after sterilization without opening the door
- Digital Display of all the operations in sequence and timer & Temperature.
- The internal chamber have two number of trays made of stainless steel
- One Sterilization Cycle with dry cycle completes with single touching of start button.
- Fitted with Silicon Gasket, Inner Aluminium & Outer M.S. duly powder coated

Front Loading Table Top Autoclave

Inner Chamber Size	Ltrs (Approx.)
225 x 310 mm	12
255 x 310 mm	16
310 x 310 mm	22

2/3, DSIDC, Community Works Centre, J-Block, Shakurpur, (Opp. Britannia Ind. Ltd.) Delhi-110034 (India)

Phone: 0091-11-27101383, 27184490, Fax: 0091-11-27108890

E-mail: info@opticstechnology.in, opticstechnology@gmail.com

Website: www.opticstechnology.in