

SINCE 1978



# OPTICS TECHNOLOGY

Manufacturer & Exporters of : LABORATORY TESTING INSTRUMENTS



An ISO 9001:2008 CERTIFIED  
An ISO 10002:2004 CERTIFIED

**CE Certified**

## GLASS BEADS STERILIZER



GLASS BEADS STERILIZER

Glass Bead sterilizers are very convenient to use over alcohol / flame for sterilizing implements (e.g. forceps, scalpels etc.,) in the laboratory.

### FEATURES:

Quick and clean alternative to conventional methods of sterilization, ensure total sterilization by destruction of anti micro organisms.

Operating temperature upto 250-265°C.

Fully insulated body.

The unit is compact and can be placed on any Laminar Air Flow Work tables in a clean room atmosphere.