

SINCE 1978



# OPTICS TECHNOLOGY

Manufacturer & Exporters of : LABORATORY TESTING INSTRUMENTS



An ISO 9001:2008 CERTIFIED  
An ISO 10002:2004 CERTIFIED

**CE Certified**

## Gyratory Sieve Shaker



Gyratory Sieve Shaker

FOR 8", 12" & 18" TEST SIEVES Gyratory Sieve Shaker is driven by ¼ H.P. Electric motor through a reduction gear, can carry upto SIX, 6" OR 8" /12" /18" Test Sieves with one set of LID AND Pan. In this Gyratory motion of the sieve shaker there is an upward and downward movement but is inclined from the vertical axis and the direction or inclination changes progressively in clockwise direction having frequency of about 270 vibration per minutes. Provided with or without timer as per requirement.

2/3, DSIDC, Community Works Centre, J-Block, Shakurpur, (Opp. Britannia Ind. Ltd.) Delhi-110034 (India)

Phone: 0091-11-27101383, 27184490, Fax: 0091-11-27108890

E-mail: [info@opticstechnology.in](mailto:info@opticstechnology.in), [opticstechnology@gmail.com](mailto:opticstechnology@gmail.com)

Website: [www.opticstechnology.in](http://www.opticstechnology.in)