

SINCE 1978



# OPTICS TECHNOLOGY

Manufacturer & Exporters of : LABORATORY TESTING INSTRUMENTS



An ISO 9001:2008 CERTIFIED  
An ISO 10002:2004 CERTIFIED

**CE Certified**

## Hand Sterilizer (Ozone)



Hand Sterilizer (Ozone)

Ozone Hand Sterilizer is used for sterilization of hands prior to working :

Sterilization: in 5-6 seconds only  
S.S. Body: (GMP model)  
Operation Time: By Timer  
Electric Consumption: 50 watts / 220 V  $\pm$  10%  
Sterility Level: 99.9% in case of wet hands  
Approach sensor: Infra Red