

SINCE 1978



# OPTICS TECHNOLOGY



Manufacturer & Exporters of : LABORATORY TESTING INSTRUMENTS



An ISO 9001:2008 CERTIFIED  
An ISO 10002:2004 CERTIFIED

**CE Certified**

## Heating Mantle



Heating Mantle

Most suitable for labs, pharmaceuticals, chemicals glass plants etc. useful for heating inflammable liquids electric heating net is hard knitted from glass yarn (without any joint) and for durability & stands. Body of mantle is spun without joint duly painted (one piece). Mantles are logged with special grade mineral wool. Each Mantle is fitted with '**ENERGY REGULATER**' for maintaining temperature upto 300°C and with pilot lamp complete with 3 core wire and plug to work on 220 volts. A.C.

CAPACITY (ml)	WATTAGE (w)
250	150
500	200
1000	300
2000	400
3000	650
5000	850

2/3, DSIDC, Community Works Centre, J-Block, Shakurpur, (Opp. Britannia Ind. Ltd.) Delhi-110034 (India)

Phone: 0091-11-27101383, 27184490, Fax: 0091-11-27108890

E-mail: [info@opticstechnology.in](mailto:info@opticstechnology.in), [opticstechnology@gmail.com](mailto:opticstechnology@gmail.com)

Website: [www.opticstechnology.in](http://www.opticstechnology.in)