

SINCE 1978



OPTICS TECHNOLOGY

Manufacturer & Exporters of : LABORATORY TESTING INSTRUMENTS



An ISO 9001:2008 CERTIFIED
An ISO 10002:2004 CERTIFIED

CE Certified

High Precision water Bath



High Precision water Bath

SPECIFICATIONS:

- Double walled, Inner chamber made of Stainless steel.
- Outer chamber made of M.S. Duly Powder Coated.
- Gap between the two walls is filled with Glass wool.
- Temp. is controlled by DIGITAL ELECTRONIC TEMP. INDICATOR CUM CONTROLLER from 5°C above room Temp. to 95°C, Accuracy $\pm 0.5^\circ\text{C}$.
- Stirrer is provided for continuous stirring
- Front glass window on one side for visibility of the sample.
- Workable on 220 Volts 50 Hz. AC mains.

Size: 280 x 360 x 380 mm - 22 Ltrs

2/3, DSIDC, Community Works Centre, J-Block, Shakurpur, (Opp. Britannia Ind. Ltd.) Delhi-110034 (India)

Phone: 0091-11-27101383, 27184490, Fax: 0091-11-27108890

E-mail: info@opticstechnology.in, opticstechnology@gmail.com

Website: www.opticstechnology.in