

SINCE 1978



# OPTICS TECHNOLOGY

Manufacturer & Exporters of : LABORATORY TESTING INSTRUMENTS



An ISO 9001:2008 CERTIFIED  
An ISO 10002:2004 CERTIFIED

**CE Certified**

## Incinerator Portable



Incinerator Portable

Incinerator is precisely designed for quick and safe disposal of waste materials specially for pharmacy colleges, small R&D labs, small nursing homes, veterinary colleges without causing any harm to the environment in a very eco-friendly process.

### CONSTRUCTION:

Body is made of mild steel duly powder coated. Heating element is made of Kanthal wire. Ducting unit is also provided with this instrument, with water circulation in the lower chamber and burning portion in the upper chamber.

Rating of heater

Lower chamber 12 KW Temperature range 850 - 900°C Upper chamber 4.5 KW, Temperature range 1050°C.

Available in different models as per specification below.

CAPACITY (ml)	Chamber size
1 Kg.	12" x 12" x 12"
3 Kg.	12" x 12" x 18"
5 kg.;	12" x 12" x 24"

2/3, DSIDC, Community Works Centre, J-Block, Shakurpur, (Opp. Britannia Ind. Ltd.) Delhi-110034 (India)

Phone: 0091-11-27101383, 27184490, Fax: 0091-11-27108890

E-mail: [info@opticstechnology.in](mailto:info@opticstechnology.in), [opticstechnology@gmail.com](mailto:opticstechnology@gmail.com)

Website: [www.opticstechnology.in](http://www.opticstechnology.in)