

SINCE 1978



OPTICS TECHNOLOGY

Manufacturer & Exporters of : LABORATORY TESTING INSTRUMENTS



An ISO 9001:2008 CERTIFIED
An ISO 10002:2004 CERTIFIED

CE Certified

Microprocessor based Karl Fischer Titrator



Features:

- Auto Computation of Moisture in ppm, % & mg H₂O.
- Storage of Date, Time, Batch No. and Final Results
- 8 Digit Digital Display.
- Motorized Dispensing
- External Printer Attachment Facility (Parallel port Dot matrix)
- Microprocessor based technology

Microprocessor based Karl Fischer Titrator

Specifications:

Measuring Range 10 ppm to 100% of moisture

Moisture Content

Calculation (Automatic) % age : (0.00 to 100%)

ppm : (10 to 999999)

mg H₂O : (0.00 to 999.99)

End Point Alarm Visual LED and Audible Alarm

Wait Time For End Point Alarm Adjustable (1 to 59 second)

Keyboard 24 keys soft touch, membrane type

Sample Storage Upto 20 latest samples

Display 8 digit bright red, seven segment LED display

Clock In-built clock with date and time display

Stirrer In- built magnetic stirrer with adjustable speed

2/3, DSIDC, Community Works Centre, J-Block, Shakurpur, (Opp. Britannia Ind. Ltd.) Delhi-110034 (India)

Phone: 0091-11-27101383, 27184490, Fax: 0091-11-27108890

E-mail: info@opticstechnology.in, opticstechnology@gmail.com

Website: www.opticstechnology.in