

SINCE 1978



# OPTICS TECHNOLOGY

Manufacturer & Exporters of : LABORATORY TESTING INSTRUMENTS



An ISO 9001:2008 CERTIFIED  
An ISO 10002:2004 CERTIFIED

**CE Certified**

## Lab Centrifuge



Lab Centrifuge

These Centrifuges are ideal for growing routine application in biochemical and clinical labs and for determination of insolubility index for Dairy Milk Products. Unit has digital speed indicator and 0-15 minute digital countdown timer.

### Salient Features:

- Stepless speed regulator.
- Safety Lid interlock to prevent cover opening during centrifugation.
- Digital speed meter and 0-15 minutes digital countdown timer.
- Press Switch for quick acceleration to full speed.

### Technical details:

Max speed 5000 R.P.M.

**With Rotor Head 4 X 50 ml OR 8 X 15 ml**

OTHER MODELS ARE AVAILABLE ON REQUEST

2/3, DSIDC, Community Works Centre, J-Block, Shakurpur, (Opp. Britannia Ind. Ltd.) Delhi-110034 (India)

Phone: 0091-11-27101383, 27184490, Fax: 0091-11-27108890

E-mail: [info@opticstechnology.in](mailto:info@opticstechnology.in), [opticstechnology@gmail.com](mailto:opticstechnology@gmail.com)

Website: [www.opticstechnology.in](http://www.opticstechnology.in)