

SINCE 1978



OPTICS TECHNOLOGY

Manufacturer & Exporters of : LABORATORY TESTING INSTRUMENTS



An ISO 9001:2008 CERTIFIED
An ISO 10002:2004 CERTIFIED

CE Certified

Digital Melting Point Apparatus



Digital Melting Point Apparatus

This is a compact instrument designed for accurate determinations of Melting Point.

This instrument consists of a heater, having three small holes for holding the capillaries, one big hole for sensor.

It is provided with digital system to monitor the temperature digitally using a highly sensitive PT-100 sensor.

The range of the apparatus is upto 275° C with readability of 1° by digital display.

2/3, DSIDC, Community Works Centre, J-Block, Shakurpur, (Opp. Britannia Ind. Ltd.) Delhi-110034 (India)

Phone: 0091-11-27101383, 27184490, Fax: 0091-11-27108890

E-mail: info@opticstechnology.in, opticstechnology@gmail.com

Website: www.opticstechnology.in