

SINCE 1978



OPTICS TECHNOLOGY

Manufacturer & Exporters of : LABORATORY TESTING INSTRUMENTS



An ISO 9001:2008 CERTIFIED
An ISO 10002:2004 CERTIFIED

CE Certified

Microprocessor based Cond./TDS meter



0 to 200 ms/cm auto ranging Auto Calibration Data Storage Facility 16 X 2 LCD Display.

Microprocessor based Cond./TDS meter

PARAMETER	CONDUCTIVIITY	TDS	TEMPERATURE
Range Selection	Automatic	Automatic	
Range	0-20, 0-200uS, 0-2, 0- 20, 0-200mS	0-20, 0-200 ppm, 0-2, 0-20, 0-200 ppt0-	100°C
Accuracy	±0.5% FS ±1 Digit	±1% FS ±1 Digit	
Resolution	0.01 uS	0.1°C	
Cell constant	Adjustable (0.05-10.5)	0.01 ppm0.	1°C
Temperature compensation	Auto/Manual (0-50°C)	Auto/Manual	
Temperature Coefficient	Adjustable (0.1-4.0%)	Adjustable (0.1-4.0%)	
Readout	16x2 Line Alphanumeric LCD display		
Storage	Up to 90 Samples		
Key board	8 soft touch membrane type keys		
Printer Interface	Dot-Matrix Printer with Centronics Interface		
Power	230V±10%AC,50Hz		

2/3, DSIDC, Community Works Centre, J-Block, Shakurpur, (Opp. Britannia Ind. Ltd.) Delhi-110034 (India)

Phone: 0091-11-27101383, 27184490, Fax: 0091-11-27108890

E-mail: info@opticstechnology.in, opticstechnology@gmail.com

Website: www.opticstechnology.in