

SINCE 1978



OPTICS TECHNOLOGY

Manufacturer & Exporters of : LABORATORY TESTING INSTRUMENTS



An ISO 9001:2008 CERTIFIED
An ISO 10002:2004 CERTIFIED

CE Certified

Microprocessor based Karl Fischer Titrator



Features:

- Auto Computation of Moisture in ppm, % & mg H₂O.
- Storage of Date, Time, Batch No. and Final Results
- 8 Digit Digital Display.
- Motorized Dispensing
- External Printer Attachment Facility (Parallel port Dot matrix)
- Microprocessor based technology

Microprocessor based Karl Fischer Titrator

Specifications:

| | |
|-------------------------------|--|
| Measuring Range | 10 ppm to 100% of moisture |
| Moisture Content | |
| Calculation (Automatic) | % age : (0.00 to 100%) |
| | ppm : (10 to 999999) |
| | mg H ₂ O : (0.00 to 999.99) |
| End Point Alarm | Visual LED and Audible Alarm |
| Wait Time For End Point Alarm | Adjustable (1 to 59 second) |
| Keyboard | 24 keys soft touch, membrane type |
| Sample Storage | Upto 20 latest samples |
| Display | 8 digit bright red, seven segment LED display |
| Clock | In-built clock with date and time display |
| Stirrer | In- built magnetic stirrer with adjustable speed |

2/3, DSIDC, Community Works Centre, J-Block, Shakurpur, (Opp. Britannia Ind. Ltd.) Delhi-110034 (India)

Phone: 0091-11-27101383, 27184490, Fax: 0091-11-27108890

E-mail: info@opticstechnology.in, opticstechnology@gmail.com

Website: www.opticstechnology.in