

SINCE 1978



OPTICS TECHNOLOGY

Manufacturer & Exporters of : LABORATORY TESTING INSTRUMENTS



An ISO 9001:2008 CERTIFIED
An ISO 10002:2004 CERTIFIED

CE Certified

Ozone Hand Sterilizer



Ozone Hand Sterilizer

Ozone Hand Sterilizer is used for sterilization of hands prior to working :

Sterilization: in 5-6 seconds only
S.S. Body: (GMP model)
Operation Time: By Timer
Electric Consumption: 50 watts / 220 V \pm 10%
Sterility Level: 99.9% in case of wet hands
Approach sensor: Infra Red

2/3, DSIDC, Community Works Centre, J-Block, Shakurpur, (Opp. Britannia Ind. Ltd.) Delhi-110034 (India)

Phone: 0091-11-27101383, 27184490, Fax: 0091-11-27108890

E-mail: info@opticstechnology.in, opticstechnology@gmail.com

Website: www.opticstechnology.in