

SINCE 1978



OPTICS TECHNOLOGY

Manufacturer & Exporters of : LABORATORY TESTING INSTRUMENTS



An ISO 9001:2008 CERTIFIED
An ISO 10002:2004 CERTIFIED

CE Certified

Penskey Martin Flash Point Apparatus



Penskey Martin Flash Point Apparatus

Useful for determination of close sup flash points above 49° C. of fuels, oils, other viscous materials and suspension of solid.

It consists of a brass test jet with handle removable, cover with spring operated rotating shutter having pilot test jet flame device, stirrer with flexible shaft.

The assembly is rested on a air bath which is covered with Dome Shape metal top.

The unit comes with a round shaped heater with a suitable energy regulator and operates on 230 V Mains AC supply.

OPTIONAL

Electrically Driver Stirrer with a FHP Motor having flexible shaft.

Spare Oil Cup Thermometer (-70°C to 110°C)

Water Bath Thermometer (-90°C to 370°C)

2/3, DSIDC, Community Works Centre, J-Block, Shakurpur, (Opp. Britannia Ind. Ltd.) Delhi-110034 (India)

Phone: 0091-11-27101383, 27184490, Fax: 0091-11-27108890

E-mail: info@opticstechnology.in, opticstechnology@gmail.com

Website: www.opticstechnology.in