

SINCE 1978



OPTICS TECHNOLOGY

Manufacturer & Exporters of : LABORATORY TESTING INSTRUMENTS



An ISO 9001:2008 CERTIFIED
An ISO 10002:2004 CERTIFIED

CE Certified

Polishing Machine



Polishing Machine

Polishing (for Metallurgical Samples) machines are used for preparing metallographic specimens & petrographic from pregrinding to smooth surface finish, for microstructure examination.

The machine is basically of motor driven disc type to provide stepless variable speed from 400 to 1100 RPM approx. Provided with a suitable H.P. single phase 220V A.C. Motor.

The brass topped disc of 8".

Available in single disc OR double disc model

2/3, DSIDC, Community Works Centre, J-Block, Shakurpur, (Opp. Britannia Ind. Ltd.) Delhi-110034 (India)

Phone: 0091-11-27101383, 27184490, Fax: 0091-11-27108890

E-mail: info@opticstechnology.in, opticstechnology@gmail.com

Website: www.opticstechnology.in