

SINCE 1978



# OPTICS TECHNOLOGY

Manufacturer & Exporters of : LABORATORY TESTING INSTRUMENTS



An ISO 9001:2008 CERTIFIED  
An ISO 10002:2004 CERTIFIED

**CE Certified**

## Pressure Module



Pressure Module

Positive pressure modules are designed to provide a clean area within a small enclosure or cell. These pressure module adopts the same Laminar Air Flow Principle. It generate positive pressure where the Room is totally Closed . It provides a source of Fresh Air in the Room.

### FEATURES

Wooden/M.S Duly Power Coated/Stainless Steel 304 Body.  
HEPA Filter ,Pre-Filter ,Blower (Motor Blower Assembly) , U.V. Lamp,  
Operating Switch  
High quality HEPA filter provides protection for the product (filtration efficiency of more than 99.97% at 0.3 micron.  
Ultra-quiet noise level of not more than 65 dBA.  
Static Pressure Manometer.  
Air Velocity :  $90 \pm 20$  FPM (Average)  
Prefilter: 10 microns  
Power Supply  
220V Single Phase, 50 Hz.

Standard Size: 2' x 2' 3' x 2' 4' x 2' 6' x 2' 8' x 2'

2/3, DSIDC, Community Works Centre, J-Block, Shakurpur, (Opp. Britannia Ind. Ltd.) Delhi-110034 (India)

Phone: 0091-11-27101383, 27184490, Fax: 0091-11-27108890

E-mail: [info@opticstechnology.in](mailto:info@opticstechnology.in), [opticstechnology@gmail.com](mailto:opticstechnology@gmail.com)

Website: [www.opticstechnology.in](http://www.opticstechnology.in)