

SINCE 1978



OPTICS TECHNOLOGY

Manufacturer & Exporters of : LABORATORY TESTING INSTRUMENTS



An ISO 9001:2008 CERTIFIED
An ISO 10002:2004 CERTIFIED

CE Certified

Redwood Viscometer



Redwood Viscometer

These viscometers are designed for viscosity tests of petroleum Products. They conform to requirements of IP 70.

Available models:

Redwood suitable for liquids having redwood flow between 20 to 2000 seconds.

Redwood suitable for liquids having redwood flow exceed 2000 seconds.

2/3, DSIDC, Community Works Centre, J-Block, Shakurpur, (Opp. Britannia Ind. Ltd.) Delhi-110034 (India)

Phone: 0091-11-27101383, 27184490, Fax: 0091-11-27108890

E-mail: info@opticstechnology.in, opticstechnology@gmail.com

Website: www.opticstechnology.in