

SINCE 1978



OPTICS TECHNOLOGY

Manufacturer & Exporters of : LABORATORY TESTING INSTRUMENTS



An ISO 9001:2008 CERTIFIED
An ISO 10002:2004 CERTIFIED

CE Certified

Rotary Flask Shaker



Rotary Flask Shaker

Rotator shaker is used in blood banks and hematology laboratories for blood typing, cross matching and mixing of donor blood with anti-coagulant, V.D.R.L. Tests and other test.

The complete body is mounted on heavy mild steel angle iron frame making it rigid and sturdy for smooth operations. Rotates specimens in horizontal plane in 3/4" circles, constructed of mild steel finished with stoved enamel point Turned up edges and rubber sheeting on platform to prevent slipping of specimen containers Platform rigidly supported by steel legs for smooth and silent rotation, mounted on four rubber suction feet to prevent creeping

- Available in mild steel platform fitted with rubber disk holders OR stainless steel lotus type clamps to hold flasks.
- Fitted with heavy duty motor having variable speed.
- Speed can be controlled through a continuously variable speed regulator knob.
- Available in different configurations of platform size, capacity and number of flasks

DIFFERENT PLATFORM SIZES AVAILABLE:-

Size:	:	
12" X 12"	:	(9 Flasks of 250 ml)
18" X 18"	:	(16 Flasks of 250 ml)
24" X 24"	:	(25 Flasks of 250 ml)
36" X 36"	:	(36 Flasks of 250 ml)
48" x 48"	:	(Capacity of flasks and qty in Nos as per requirement)
60" x 60"	:	(Capacity of flasks and qty in Nos as per requirement)
OR ANY OTHER NON - STANDARD SIZE		
OPTIONAL :		
Digital RPM Indicator		
Digital Timer		

2/3, DSIDC, Community Works Centre, J-Block, Shakurpur, (Opp. Britannia Ind. Ltd.) Delhi-110034 (India)

Phone: 0091-11-27101383, 27184490, Fax: 0091-11-27108890

E-mail: info@opticstechnology.in, opticstechnology@gmail.com

Website: www.opticstechnology.in