

SINCE 1978



OPTICS TECHNOLOGY

Manufacturer & Exporters of : LABORATORY TESTING INSTRUMENTS



An ISO 9001:2008 CERTIFIED
An ISO 10002:2004 CERTIFIED

CE Certified

Seed Germinator



Seed Germinator

Most suitable for conducting various experiments on variety of seeds under different conditions of temperature and humidity. Atmospheric conditions are created within the chamber.

- Temperature from 5°C to 50° C \pm 1°C Humidity range Ambient to 95% \pm 3%.
- Temp. & Humidity is controlled by Digital Temp. cum Humidity Controller.
- Each chamber has 14 trays.
- The Instrument has stainless steel water Reservoir at back for humidity.
- Double walled inside made of Stainless Steel and Outside made of Mild Steel with Powder Coating.
- External lighting arrangement by tube lights and Mechanical Cyclic Timer 0-24 hrs.
- Work to on 220/230V A.C.

Dimension of each chamber (Inside)
(WXDXH) 600 X 600 X 900 mm (Approx.)
24" X 24" X 36" (Approx.)

OPTIONAL:-

- A) Automatic cyclic timer (0-24 Hrs.)
- B) Dual setting digital temperature controller for setting two different temperatures for day and night.
- C) Automatic voltage stabilizer
- D) Microprocessor based Temp. & Humi. Controller



NOTE: DUAL CHAMBER MODEL & OUTER STAINLESS STEEL MODELS ARE ALSO AVAILABLE TOUCH SCREEN DISPLAY MODELS WITH DATA LOGGING OPTION ARE ALSO AVAILABLE ON REQUEST.

2/3, DSIDC, Community Works Centre, J-Block, Shakurpur, (Opp. Britannia Ind. Ltd.) Delhi-110034 (India)

Phone: 0091-11-27101383, 27184490, Fax: 0091-11-27108890

E-mail: info@optics technology.in, optics technology@gmail.com

Website: www.optics technology.in