

SINCE 1978



OPTICS TECHNOLOGY

Manufacturer & Exporters of : LABORATORY TESTING INSTRUMENTS



An ISO 9001:2008 CERTIFIED
An ISO 10002:2004 CERTIFIED

CE Certified

Slide Warming Table



Slide Warming Table

This is a heating block unit with provision for different slides at a time. Temperature Range being Amb. +5 to 250oC within $\pm 1\%$ accuracy with utmost uniformity of temperature. These are fitted with microprocessor based PID Temperature controllers

Available in models of

M.S. Body with Cast Iron Top

OR

M.S. Body with S.S. Top

OR

S.S. Body with S.S. Top

Size (In Inches)

10" X 12" X 6" (250 x 300 x 150 mm)

10" x 16" X 6" (250 x 400 x 150 mm)

12" X 18" X 6" (300 x 450 x 150 mm)

18" X 24" X 6" (450 x 600 x 150 mm)

2/3, DSIDC, Community Works Centre, J-Block, Shakurpur, (Opp. Britannia Ind. Ltd.) Delhi-110034 (India)

Phone: 0091-11-27101383, 27184490, Fax: 0091-11-27108890

E-mail: info@opticstechnology.in, opticstechnology@gmail.com

Website: www.opticstechnology.in