

SINCE 1978



# OPTICS TECHNOLOGY

Manufacturer & Exporters of : LABORATORY TESTING INSTRUMENTS



An ISO 9001:2008 CERTIFIED  
An ISO 10002:2004 CERTIFIED

**CE Certified**

## Static Pass Box



Static Pass Box

The pass box is an equipment normally installed on the wall between the clean room and any other room. It is used for transfer the material from any room to clean room or viceversa, without contaminating the air of clean room and without having to open room door. This is achieved with the help of interlocking type doors of pass box.

### Construction :

Body is made of Wooden/ M.S. Duly Powder coated/S.S. 304.

Each hinged door has wide glass viewing window, through which material kept in the pass box can be seen.

These doors are mechanical interlocked type, so only one door can be opened at a time.

U.V. Light & fluorescent light is provided in the chamber.

### Electronic Inter-Locking system

WORKING AREA SIZE :

W x D x H

2' x 2' x 2'

Other Sizes on request

2/3, DSIDC, Community Works Centre, J-Block, Shakurpur, (Opp. Britannia Ind. Ltd.) Delhi-110034 (India)

Phone: 0091-11-27101383, 27184490, Fax: 0091-11-27108890

E-mail: [info@opticstechnology.in](mailto:info@opticstechnology.in), [opticstechnology@gmail.com](mailto:opticstechnology@gmail.com)

Website: [www.opticstechnology.in](http://www.opticstechnology.in)