

SINCE 1978



OPTICS TECHNOLOGY

Manufacturer & Exporters of : LABORATORY TESTING INSTRUMENTS



An ISO 9001:2008 CERTIFIED
An ISO 10002:2004 CERTIFIED

CE Certified

Stirrer (Lab)



Stirrer (Lab)

The Motor stirrer having maximum speed upto 3000 RPM unit is clamped to a chrome plated rod fitted on a heavy cast iron base. The height and angle of the stirrer can be adjusted according to the requirement. Variable speed motor with speed regulator is fitted with S.S. stirrer rod to give smooth and noiseless operation. Available with or without Drill Chuck.

Available in following models:

- For 1 & 2 Ltr. (1/20 H.P. Motor)
- For upto 5 Ltr. (1/8 H.P. motor)
- For upto 10 Ltr. (1/4 H.P. Motor)

OPTIONAL

Digital RPM Indicator

2/3, DSIDC, Community Works Centre, J-Block, Shakurpur, (Opp. Britannia Ind. Ltd.) Delhi-110034 (India)

Phone: 0091-11-27101383, 27184490, Fax: 0091-11-27108890

E-mail: info@opticstechnology.in, opticstechnology@gmail.com

Website: www.opticstechnology.in