

SINCE 1978



OPTICS TECHNOLOGY

Manufacturer & Exporters of : LABORATORY TESTING INSTRUMENTS



An ISO 9001:2008 CERTIFIED
An ISO 10002:2004 CERTIFIED

CE Certified

Tablet Disintegration Test Apparatus



Tablet Disintegration Test Apparatus

The instrument is designed for accurate estimation of disintegration time of tablets/capsules. It is assembled in painted sheet metal cabinet with chemical resistant epoxy coating, fitted with superior quality constant speed sturdy motor. The guide discs and strokes are conforming to Pharmacopoeia Standards.

SPECIFICATIONS:

Displacement : 55mm
RPM : 30
Drive : 1/8 HP
Heater : 125 watt.
Power : AC/230 volt Single Phase.

Available in SINGLE AND DOUBLE BASKET MODELS WITH OUTER BODY OF CONSTRUCTION BEING M.S.DULY POWDER COATED BODY OR S.S. 304.

Optional :

- 1) Digital Time Totalizer
- 2) Water bath System Instead of Hot Plate.
- 3) Digital Temperature Controller in Water bath System only.



2/3, DSIDC, Community Works Centre, J-Block, Shakurpur, (Opp. Britannia Ind. Ltd.) Delhi-110034 (India)

Phone: 0091-11-27101383, 27184490, Fax: 0091-11-27108890

E-mail: info@opticstechnology.in, opticstechnology@gmail.com

Website: www.opticstechnology.in