

SINCE 1978



OPTICS TECHNOLOGY

Manufacturer & Exporters of : LABORATORY TESTING INSTRUMENTS



An ISO 9001:2008 CERTIFIED
An ISO 10002:2004 CERTIFIED

CE Certified

Tablet Dissolution Rate Test Apparatus



Tablet Dissolution Rate Test Apparatus

Dissolution Rate Test Apparatus with water bath (Single Stage)

The apparatus is suitable for testing of Dissolution rate of tablet and capsule, it consists of main transparent acrylic tank housing. Thermostatically/ Digitally controlled immersion heater, drain cock and transparent hemispherical unbreakable vessel. The base unit is made of chemical resistant powder coated sturdy frame housing, Control panel which includes Electronic regulator and thermostat. An electronic motor is mounted on chrome plated rod, operates Paddle/Basket.

SPECIFICATIONS:

MOTOR : Electric Motor "
HEATER : 500 Watts.
TEMP. : $37.0 \pm 0.50^{\circ}\text{C}$
POWER : AC/220 Volts / Single Phase

Optional :

1. Digital RPM Indicator
2. Digital Timer
3. Digital temperature indicator

2/3, DSIDC, Community Works Centre, J-Block, Shakurpur, (Opp. Britannia Ind. Ltd.) Delhi-110034 (India)

Phone: 0091-11-27101383, 27184490, Fax: 0091-11-27108890

E-mail: info@opticstechnology.in, opticstechnology@gmail.com

Website: www.opticstechnology.in