

SINCE 1978



OPTICS TECHNOLOGY

Manufacturer & Exporters of : LABORATORY TESTING INSTRUMENTS



An ISO 9001:2008 CERTIFIED
An ISO 10002:2004 CERTIFIED

CE Certified

Test Tube Warmer (Specially for I.V.F. Labs)



Test Tube Warmer (Specially for I.V.F. Labs)

Digital units include built-in temperature sensing probe Compact design saves valuable bench space. Powder-coated steel body construction ensures durability. Digital units utilize a Digital Temperature Controller and easy-to-read LED display for accurate and reproducible temperature control.

Temperature Range: ambient +5°C to 50°C.

S.S. body with anodized aluminium blocks for heating test tubes of size 12 x 13 ml tubes OR 12 x 6 ml tubes

2/3, DSIDC, Community Works Centre, J-Block, Shakurpur, (Opp. Britannia Ind. Ltd.) Delhi-110034 (India)

Phone: 0091-11-27101383, 27184490, Fax: 0091-11-27108890

E-mail: info@opticstechnology.in, opticstechnology@gmail.com

Website: www.opticstechnology.in