

SINCE 1978



# OPTICS TECHNOLOGY



Manufacturer & Exporters of : LABORATORY TESTING INSTRUMENTS



An ISO 9001:2008 CERTIFIED  
An ISO 10002:2004 CERTIFIED

**CE Certified**

## Tintometer



- \* Visual measurement of color matching to determine material color.
- \* Design in conformity with internationally accepted color scale of Lovibond.
- \* Simple structure and easy to operate.
- \* With different size cuvettes and better observe the visual color effect of various materials.

### Specifications:

Tintometer

Measurement range	Red	Yellow	Blue	Neutral
in Lovibond units	R0.1 – R79.9	Y0.1 – Y79.9	B0.1 – B49.9	N0.1 – N3.9
Minimum Readout	0.1 Lovibond Unit			
Whiteboard reflectivity	>80%			
Magnification	1.9X			
Cuvette Size	10 mm x 20 mm x 40 mm, 25.4 mm x 20 mm x 40 mm 133.4 mm x 20 mm x 40 mm			
Power Supply	230 V ± 22V	50Hz		

2/3, DSIDC, Community Works Centre, J-Block, Shakurpur, (Opp. Britannia Ind. Ltd.) Delhi-110034 (India)

Phone: 0091-11-27101383, 27184490, Fax: 0091-11-27108890

E-mail: info@opticstechnology.in, opticstechnology@gmail.com

Website: www.opticstechnology.in