

SINCE 1978



OPTICS TECHNOLOGY

Manufacturer & Exporters of : LABORATORY TESTING INSTRUMENTS



An ISO 9001:2008 CERTIFIED
An ISO 10002:2004 CERTIFIED

CE Certified

Tissue Floatation Bath



Tissue Floatation Bath

Temperature Range ambient +5deg C to 80 deg C.

- Interior of the bath made of black Anodized Aluminum.
- Digital temperature Controller-cum-Indicator
- LED display for Set Value & Process Value

Tissue floatation bath size: minimum 200 mm X 100 mm

2/3, DSIDC, Community Works Centre, J-Block, Shakurpur, (Opp. Britannia Ind. Ltd.) Delhi-110034 (India)

Phone: 0091-11-27101383, 27184490, Fax: 0091-11-27108890

E-mail: info@opticstechnology.in, opticstechnology@gmail.com

Website: www.opticstechnology.in