

SINCE 1978



OPTICS TECHNOLOGY

Manufacturer & Exporters of : LABORATORY TESTING INSTRUMENTS



An ISO 9001:2008 CERTIFIED
An ISO 10002:2004 CERTIFIED

CE Certified

Weighing Machine Electronic



- External calibration / Motorized Internal Calibration with built in weight
- A built in under hook for below weighing
- Easy to read Backlight LCD Display
- Density measurement for solid & liquid (Optional)
- Standard RS 232 C Interface (Optional)
- Auto - Timing Calibration
- Percentage (%) Weighing
- High Resolution, Quick Weighing,
- Based on E.M.F.C. Technology
- Tare range upto max. Capacity
- Multiple Weighing Units
- Level Indicator, Adjustable fit
- Counting Weighing mode
- Units Conversion and data output
- Overload Protection

Weighing Machine Electronic

CAPACITY	:	220 Gms
READABILITY	:	0.1 mg.
REPEATABILITY (±)	:	0.1 mg.
LINEARITY	:	0.2 mg.
PAN SIZE (mm/inch)	:	90 mm
RESPONSE TIME	:	3-4 sec.
CALIBRATION	:	External / Motorized Internal
TARE RANGE	:	Full
OPERATING TEMP.	:	15°C to 35°C
POWER SUPPLY	:	A/C adapter 230 V

NOTE: ELECTRONIC BALANCES IN OTHER CAPACITIES WITH DIFFERENT READABILITY ARE ALSO AVAILABLE ON REQUEST.

2/3, DSIDC, Community Works Centre, J-Block, Shakurpur, (Opp. Britannia Ind. Ltd.) Delhi-110034 (India)

Phone: 0091-11-27101383, 27184490, Fax: 0091-11-27108890

E-mail: info@opticstechnology.in, opticstechnology@gmail.com

Website: www.opticstechnology.in