

SINCE 1978



# OPTICS TECHNOLOGY

Manufacturer & Exporters of : LABORATORY TESTING INSTRUMENTS



An ISO 9001:2008 CERTIFIED  
An ISO 10002:2004 CERTIFIED

**CE Certified**

## Wrist Action Shaker



Wrist Action Shaker

The shaking machine fitted with D.C. MOTOR to work on 220 VOLT A/C 50 CYCLES to work for longer period smoothly. Flasks mounted spindle to oscillate through a small amplitude sturdy eccentric mechanism provides effective mechanical shaking . The shaking machine constructed to hold 100 ml to 1000 ml. Conical flasks. The shaking speed action is controlled by ELECTRONIC regulator . 100 R.P.M TO 800 R.P.M approx

### HOLDING CAPACITY

- (A) To hold 4 Flasks Cap. 100 ml. to 1000 ml. by side      clamps
- (B) To hold 8 Flasks Cap. 100 ml. to 1000 ml. by side      clamps

2/3, DSIDC, Community Works Centre, J-Block, Shakurpur, (Opp. Britannia Ind. Ltd.) Delhi-110034 (India)

Phone: 0091-11-27101383, 27184490, Fax: 0091-11-27108890

E-mail: [info@opticstechnology.in](mailto:info@opticstechnology.in), [opticstechnology@gmail.com](mailto:opticstechnology@gmail.com)

Website: [www.opticstechnology.in](http://www.opticstechnology.in)