



# OPTICS TECHNOLOGY



SINCE 1978

Manufacturer & Exporters of : LAB TESTING INSTRUMENTS



An ISO 9001:2015 CERTIFIED  
An ISO 10002:2018 CERTIFIED  
An ISO 13485:2016 CERTIFIED

IEC 61010-1:2010  
BS EN 12469:2000  
ISO 29463-5:2022

## LAB STIRRER WITH VARIABLE SPEED & S.S. SHAFT

CE



OT-LS-53 MS



OT-LS-53 SS

GMP Model

The motor stirrer having maximum speed upto 3000 RPM unit is clamped to a Chrome plated rod fitted on a heavy cast iron base. The height and angle of the stirrer can be adjusted according to the requirement. Variable speed motor with speed regulator is fitted with S.S. Stirrer rod to give smooth and Noiseless operation. **Available with Drill-Chuck.**

Model No.	Available in following models
OT-LS-53/MS/SS-A	: For 1Ltr. (1/20 H.P.Motor)
OT-LS-53/MS/SS-B	: For 2 Ltr. (1/20 H.P.Motor)
OT-LS-53/MS/SS-C	: For 5 Ltr. (1/8 H.P.Motor )
OT-LS-53/MS/SS-D	: For upto 10 Ltr. (1/4 H.P.Motor)
OT-LS-53/MS/SS-E	: For upto 20 Ltr. <b>OR</b> any other capacity on request

### Optional :

- 1) Digital RPM Indicator
- 2) Digital Timer

2/3, DSIDC, (Behind Fun Cinema) NSP, Pitampura, Delhi-110034 (INDIA)

☎ 9811238150, 9811238151, 9811236960

E-mail: [info@opticstechnology.in](mailto:info@opticstechnology.in),

[opticstechnology@gmail.com](mailto:opticstechnology@gmail.com),

Website : [www.opticstechnology.in](http://www.opticstechnology.in)