

SINCE 1978



OPTICS TECHNOLOGY

Manufacturer & Exporters of : LABORATORY TESTING INSTRUMENTS



An ISO 9001:2008 CERTIFIED
An ISO 10002:2004 CERTIFIED

CE Certified

Ball Mill



Ball Mill

For mixing & grinding, electrically operated having capacity of few gms to 1 Kg & 2 Kg. Fitted with geared Motor with max. speed of 80 rpm. The jar is made of SS and having steel balls of different sizes to be put in jar for mixing & grinding pesticides powder. Jar is interchangeable and easily removable, light in weight and compact unit working on 220 V AC with cap. 1 KG. Variable speed regulator is provided

Capacity 1 kg
Capacity 2 kg
Capacity 5 kg
Capacity 10 kg
Capacity 20 Kg

2/3, DSIDC, Community Works Centre, J-Block, Shakurpur, (Opp. Britannia Ind. Ltd.) Delhi-110034 (India)

Phone: 0091-11-27101383, 27184490, Fax: 0091-11-27108890

E-mail: info@opticstechnology.in, opticstechnology@gmail.com

Website: www.opticstechnology.in