

SINCE 1978



# OPTICS TECHNOLOGY

Manufacturer & Exporters of : LABORATORY TESTING INSTRUMENTS



An ISO 9001:2008 CERTIFIED  
An ISO 10002:2004 CERTIFIED

**CE Certified**

## Microbiological Blender



Microbiological Blender

(With Autoclavable Jar)

### TECHNICAL SPECIFICATIONS

RPM	14000- 18000 RPM Approx.
Watts	450 Watts (2 Amp.)
Timer	Selective Timer Cycle from 0.1 to 5 minutes for automatic blending or manual control setting for continuous running.
Base	Aluminium Cast - Baked Enamel Finish
Container	Stainless Steel Jar with S.S. Blades & S.S. Cap.
Capacity	(a) 250ml (b) 500 ml. (c) 1000ml.
Optional	50-100 ml containers of Stainless Steel.

2/3, DSIDC, Community Works Centre, J-Block, Shakurpur, (Opp. Britannia Ind. Ltd.) Delhi-110034 (India)

Phone: 0091-11-27101383, 27184490, Fax: 0091-11-27108890

E-mail: [info@opticstechnology.in](mailto:info@opticstechnology.in), [opticstechnology@gmail.com](mailto:opticstechnology@gmail.com)

Website: [www.opticstechnology.in](http://www.opticstechnology.in)