

SINCE 1978



OPTICS TECHNOLOGY

Manufacturer & Exporters of : LABORATORY TESTING INSTRUMENTS



An ISO 9001:2008 CERTIFIED
An ISO 10002:2004 CERTIFIED

CE Certified

Clean Air Tent



Clean Air Tent

These are overhead clean air modules mounted on roof or with steel frame on the ground. Units can be built with / without casters so the clean enclosure can be moved over equipment or used to isolate a controlled process. It is basically used for filling, refilling weighing and sampling/Dispensing of Raw Material, Basic Compounds, Active Compounds, works on the principle of laminar air flow with HEPA Filters & Pre-Filters. Available in Wooden / M.S. body duly powder coated / Stainless steel body.

Specifications:-

Standard Size (Working Area)	2' x 2', 3' x 3', 4' x 4', 6' x 6', 8' x 8'
Air Velocity	90 ± 20 cfm (Average)
Filtration Elements	Main Filters : HEPA Filter
Prefilter:	Washable non-woven polyester fibers
Filter Efficiency	Main Filters : 99.97 at 0.3 micron Prefilter: Arrestance 90%, at 10 micron
Noise Level	< 65 dBA
Light Intensity	800-1000 Lux
Power Supply	220V Single Phase, 50 Hz.

2/3, DSIDC, Community Works Centre, J-Block, Shakurpur, (Opp. Britannia Ind. Ltd.) Delhi-110034 (India)

Phone: 0091-11-27101383, 27184490, Fax: 0091-11-27108890

E-mail: info@opticstechnology.in, opticstechnology@gmail.com

Website: www.opticstechnology.in